





SYSTEM REPORT: Phillips Access 16

Project No. (Item)	004508
	QA
Inspection done by & Date	B.Roosdorp 29-sept-2022
System Information	
Manufacturer and model	Philips Access
Upgraded System? <ul style="list-style-type: none"> • From which one? • Year of the upgrade? 	Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>
System model no.	728116
System serial no.	310175
System YOM	2018
Number of slices	16
Air or water cooling?	Air
Gantry & Detector Information	
Gantry YOM	2018
Gantry model no.	CTV16PFA100
Gantry serial no.	302310
Detector model no.	
Tube Information	
Tube Type	CTR1735
Tube Housing model no.	4598-005-94291
Tube Housing serial no.	17471K277687
Tube Insert model no.	DU 3506
Tube Insert serial no.	277687
Tube YOM	2021 -11
Tube installation date	22-12-2021

Patient Table Information	
Table model no.	CTV16PVA200
Table serial no.	306209
Table YOM	2018
Head rest present?	Yes <input checked="" type="checkbox"/> NO <input type="checkbox"/>
Positioning pads incl.?	Yes <input checked="" type="checkbox"/> NO <input type="checkbox"/>
WORKSTATION INFORMATION (Main-Host)	
Workstation model no.	CTV16PFA300
Workstation serial no.	306300
Workstation software version	2.0
Workstation software options	BOLUS TRACKING, IDOSE4, WORKLIST, MPPS
ACCESSORIES, MANUAL HARD-COPIES AND SOFTWARE DISCS	
Software Manuals	Yes <input checked="" type="checkbox"/> NO <input type="checkbox"/> Yes <input checked="" type="checkbox"/> NO <input type="checkbox"/>
Any additional items/accessories?	1. Control box gantry 
MAINTENANCE	
Any major part replacements/service done?	1. None recent according to the user. 2. Last maintenance not done as system was not in use

COMMENTS	
Functional condition	<p>1. Normal</p> <p>Patient data & hospital information will be removed by the operator before de-install begins</p>
Cosmetic condition	No scratches.

PHANTOMS		
Model name	Model no.	Serial no.
Head	455012201471 - H	180508337
Body	No label	
Water	No label	

Q&A TESTING		
	Tests Names Performed	RESULT OF TESTS
		QA
1	Daily QA Pass	Pass